




Basic knowledge in plasma surgery


Argon plasma coagulation



Principles of plasma surgery



Schematic illustration of monopolar argon plasma coagulation



The argon gas is ignited at the electrode of the APC probe, the ionized argon plasma transfers the energy to the target tissue

ARGON PLASMA COAGULATION

Argon plasma coagulation (APC) is an electrosurgical procedure. Here, high frequency alternating current is transferred to the target tissue from the tip of the probe via the ionized argon gas. This procedure reliably arrests bleeding with an effective and metered surface coagulation and devitalizes tissue. The APC is contact-free so that the distal end of the instrument cannot adhere to the coagulated tissue and tear open the scab that has formed. A further advantage is the limited penetration depth of the APC which minimizes the risk of perforations.¹

Due to its numerous application advantages, the procedure is employed in endoscopy and open surgery.

PHYSICAL PRINCIPLES²⁻⁴

In APC, the energy is transferred to the tissue via ionized, electrically conductive argon plasma using an APC probe with monopolar technology. The thermal effects are coagulation, desiccation or devitalization of the target tissue.

In contrast to lasers, the energy between the electrode and the target tissue is transferred via an electrical field in APC and not optically. The argon plasma beam follows the path of least electrical resistance.

¹ Kähler, G F et al. Investigation of the thermal tissue effects of the argon plasma coagulation modes pulsed and precise on the porcine esophagus, ex vivo and in vivo; *Gastrointest. Endosc.*, 2009

² Eickhoff A, Repici A, Manner H, Enderle, MD. *Electrosurgical Pocket Guide for GI Interventions*; Erbe Elektromedizin GmbH


³ Eickhoff, A et al.: Prospective nonrandomized comparison of two modes of argon beamer (APC) tumor desobstruction: effectiveness of new pulsed APC versus forced APC; *Endoscopy* 2007

⁴ Zenker, M. *Argon plasma Coagulation*; GMS Krankenhhyg Interdisziplin. 2008

Tissue effects

The tissue effect of the APC is provided by the current which flows through the tissue and the resulting endogenous heating. One distinguishes between various thermal effect zones in the tissue, depending on the achieved target temperature.

1. Hyperthermia, 2. Ablation 3. Coagulation/Dessiccation, 4. Carbonization and 5. Vaporization.



	above approx.
1 Hyperthermia	40° C
2 Ablation	60° C
3 Coagulation/Dessiccation	100° C
4 Carbonization	150° C
5 Vaporization	300° C

The tissue effect spreads radially in depth.⁴

FACTORS INFLUENCING THE TISSUE EFFECT^{1,4}

The degree of the thermal effect of the APC on the tissue depends on several factors. The most important factors on coagulation depth in order of priority:


1. Application duration (especially for static application)
2. Power setting or effect stage
3. Probe distance (working distance)
4. Other factors: type of tissue, static/dynamic application

1 Kähler, G F et al. Investigation of the thermal tissue effects of the argon plasma coagulation modes pulsed and precise on the porcine esophagus, ex vivo and in vivo; *Gastrointest. Endosc.*, 2009


4 Zenker, M. Argon plasma Coagulation; *GMS Krankenhhyg Interdiszip.* 2008

Tissue effects

Duration of application



Power



APPLICATION TIME – THE MOST IMPORTANT FACTOR²⁻⁷


The longer the APC is activated, the deeper the effect on the target tissue. For this reason, we recommend starting with short activation times and to increase the duration step by step up until the desired effect is achieved. In case of longer APC application at a single point, the depth effect increases strongly and if the application duration is too long, the tissue may be carbonized and perforated.

In dynamic application, the APC probe should be moved under visual control in slow, controlled movements (brushstrokes) over the target tissue.

POWER SETTING²⁻⁷

The coagulation depth is dependent on the power setting and should be set according to the localization and indication.

Probe distance



PROBE DISTANCE²⁻⁷

The greater the probe distance, the lower the penetration depth. If the probe distance increases, a point can be reached where ignition is no longer possible.

FURTHER FACTORS: TYPE OF TISSUE²⁻⁷

The structures of biological tissue differ in their sensitivity, which has to be considered with electrosurgery and particularly with APC in the power setting and application duration.

- 2 Eickhoff A, Repici A, Manner H, Enderle, MD. *Electrosurgical Pocket Guide for GI Interventions*; Erbe Elektromedizin GmbH
- 3 Eickhoff, A et al. *Prospective nonrandomized comparison of two modes of argon beamer (APC) tumor desobstruction: effectiveness of new pulsed APC versus forced APC*; *Endoscopy* 2007
- 4 Zenker, M. *Argon plasma Coagulation*. *GMS Krankenhhyg Interdiszip*. 2008
- 5 Taghavi SA, Soleimani SM, Hosseini-Asl SM et al. *Adrenaline injection plus argon plasma coagulation versus adrenaline injection plus hemoclips for treating high-risk bleeding peptic ulcers: a prospective, randomized trial*. *Can J Gastroenterol* 2009; 23(10): 699 – 704.
- 6 Wang HM, Hsu PI, Lo GH et al. *Comparison of hemostatic efficacy for argon plasma coagulation and distilled water injection in treating high-risk bleeding ulcers*. *J Clin Gastroenterol* 2009; 43(10): 941 – 945.
- 7 Herrera S, Bordas JM, Llach J et al. *The beneficial effects of argon plasma coagulation in the management of different types of gastric vascular ectasia lesions in patients admitted for GI hemorrhage*. *Gastrointest Endosc* 2008; 68(3): 440 – 446.

APC modes

The constant voltage control of the plasma modes enables consistent quality and reproducibility of the tissue effects.²²

forcedAPC



Effective devitalization with forcedAPC

preciseAPC[®]



The preciseAPC[®] mode used in thin-walled structures

pulsedAPC[®]



The pulsedAPC[®] mode used for tissue ablation and coagulation

This mode offers effective coagulation and devitalization. The HF power is voltage-controlled up to 120 Watt and is applied as a continuous energy transfer.

Unlike forcedAPC, preciseAPC[®] operates in the lower energy range. This allows uniform coagulation effects to be precisely adjusted in the target tissue, which enables a homogenous tissue effect particularly in thin-walled structures or peristaltic movement.

This mode is based on pulsed activation (on-off). pulsedAPC[®] is variable in use, for ablating or coagulating tissue. pulsedAPC[®] is easy to meter and offers homogenous tissue effects as a result. In pulsedAPC[®], powers of 1 to 120 Watt may be set. 2 different pulse frequencies can be set.



Applications


Applications in flexible endoscopy, gastroenterology and pneumology

- ☑ Chronic bleeding⁵⁻⁸ (Fig. 1)
 - GAVE syndrome (watermelon stomach)
 - Radioproctitis
 - Angiodysplasia
- ☑ Coagulation of bleeding in the resection base following EMR⁹
- ☑ Devitalization of tumor residues following EMR¹⁰
- ☑ Immediate recanalization of exophytic stenoses¹⁵⁻¹⁸
- ☑ Coagulation of diffuse and acute bleeding in the entire gastrointestinal and bronchial tract^{8,19}
- ☑ Devitalization of stent ingrowth or overgrowth²⁰
- ☑ Trimming of stents in the gastrointestinal or bronchial tract²¹



Coagulation of a telangiectasia with APC

- 5 Taghavi SA, Soleimani SM, Hosseini-Asl SM et al. Adrenaline injection plus argon plasma coagulation versus adrenaline injection plus hemoclips for treating high-risk bleeding peptic ulcers: a prospective, randomized trial. *Can J Gastroenterol* 2009; 23(10): 699 – 704.
- 6 Wang HM, Hsu PI, Lo GH et al. Comparison of hemostatic efficacy for argon plasma coagulation and distilled water injection in treating high-risk bleeding ulcers. *J Clin Gastroenterol* 2009; 43(10): 941 – 945.
- 7 Herrera S, Bordas JM, Llach J et al. The beneficial effects of argon plasma coagulation in the management of different types of gastric vascular ectasia lesions in patients admitted for GI hemorrhage. *Gastrointest Endosc* 2008; 68(3): 440 – 446.
- 8 Manner H, Enderle MD, Pech O et al. Second-generation argon plasma coagulation: two-center experience with 600 patients. *J Gastroenterol Hepatol* 2008; 23(6): 872 – 878.
- 9 Katsinelos P, Gkagkalis S, Paroutoglou G, Chatzimavroudis G, Fasoulas K, Zavos C, Varitimadakis K, Lazaraki G, Kotronis G, Kountouras J. A prospective comparative study of blended and pure coagulation current in endoscopic mucosal resection of large sessile colorectal polyps. *Surg Laparosc Endosc Percutan Tech*. 2014 Jun; 24(3): 226 – 31.
- 10 Manner H, Rabenstein T, Pech O, Braun K, May A, Pohl J, Angelika Behrens A, Vieth M, Ell C. Ablation of residual Barrett's epithelium after endoscopic resection: a randomized long-term follow-up study of argon plasma coagulation vs. surveillance (APE study). *Endoscopy* 2014; 46(01): 6–12
- 15 Wang H, Tao M, Zhang N, Luo L, Li D, Zou H, Zhou Y, Liang S. Bronchoscopic interventions combined with percutaneous modalities for the treatment of thyroid cancers with airway invasion. *Eur Arch Otorhinolaryngol*. 2015 Feb; 272(2): 445 – 51.
- 16 Wang JW, Huang M, Zha WJ, Zhou LF, Qi X, Wang H. Flexible bronchoscopic intervention for endobronchial hamartoma. *Zhonghua Jie He Hu Xi Za Zhi*. 2013 Dec; 36(12): 963 – 7. Chinese
- 17 Seaman JC, Musani AI. Endobronchial ablative therapies. *Clin Chest Med*. 2013 Sep; 34(3): 417 – 25. doi: 10.1016/j.ccm.2013.04.006.
- 18 Sim DW, Oh JJ, Kim KS, Choi YD, Kwon YS. Pleomorphic adenoma of the trachea. *J Bronchology Interv Pulmonol*. 2014 Jul; 21(3): 230 – 3.
- 19 Reichle G. Die Argonplasma-Koagulation zur bronchoskopischen Rekanalisation und Blutstillung. *Atemw- Lungenkrkh* 2003; Jahrgang 29: 258 – 269.
- 20 Reichle G, Freitag L, Kullmann HJ, Prenzel R, Macha HN, Farin G. Argon plasma coagulation in bronchology: a new method – alternative or complementary?. *Pneumologie* 2000; 54: 508 – 516.
- 21 Chen YK, Jakribettuu V, Springer EW, Shah RJ, Penberthy J, Nash SR. Safety and efficacy of argon plasma coagulation trimming of malpositioned and migrated biliary metal stents: a controlled study in the porcine model. *Am J Gastroenterol* 2006; 101: 2025 – 2030



Extensive homogenous coagulation with APC



Low smoke formation in argon-assisted cutting



Diffuse surface coagulation with APC



Extensive homogenous coagulation with APC




Good view of the operative site with argon-assisted incision

Applications in gynecology, urology and general surgery

- ☑ Coagulation of extensive bleeding in breast surgery (Fig. 1)
- ☑ Argon-assisted resection of mammary layers and mammary gland tissue (Fig. 2)
- ☑ Coagulation of a resection base in partial nephrectomy (Fig. 3)
- ☑ Coagulation of the liver bed base with the APC function of the APCapplicator (Fig. 4)
- ☑ The liver capsule is opened with an argon-assisted incision. (Fig. 5)

Erbe products* for plasma surgery

The Erbe equipment for plasma surgery consists of the workstation with VIO® 3 and APC 3 as well as open surgical, laparoscopic and endoscopic probes and applicators. The workstation supports the instruments and applications with the forcedAPC, preciseAPC® and pulsedAPC® modes. Virtually all indications can be treated with these modes, ranging from selective flat coagulation of minor bleeding to the devitalization of extensive lesions.



FiAPC® probes in various lengths and designs.



The APC applicator and its functions covers numerous procedures in the surgical specialties. The instrument is available in open surgical and laparoscopic design.

Beam forms



The FiAPC® probes have axial, lateral and circular beam forms

The electrosurgery workstation with the electrosurgery unit VIO® 3 and APC 3 on a Cart. The APC 3 is operated via the display of the VIO® 3.



THE DIRECT LINK TO THE WEBCAST
"THE FUNDAMENTALS OF ELECTROSURGERY IN ENDOSCOPY"
BY PD DR. AXEL EICKHOFF

Important information

We have prepared this document with care. Nonetheless, we cannot completely rule out errors in this document.

The information, recommendations and other data („Information“) contained in this document reflect our state of knowledge and the state of science and technology at the time of preparing the document. The information is of a general nature, non-binding and serves solely for general information purposes and does not represent instructions for use or notes on application.

The information and recommendations contained in this document do not constitute any legal obligations by Erbe Elektromedizin GmbH as well as their associated companies („Erbe“) or any other claims against Erbe. The information does not represent a guarantee or other quality statement, these require the express contractual arrangement with Erbe in individual cases. Erbe shall not be liable for any type of damage resulting from following information given in this document, regardless of the legal reason for liability.

Every user of an Erbe product is responsible for checking the respective Erbe product for its properties as well as the suitability for the intended type of application or intended purpose in advance. The suitable type of application of the respective Erbe product is given by the user manual and the notes on use for the corresponding Erbe product. The user is obliged to check whether the existing user manual and the notes on use correspond with the status for the specific Erbe product. The devices may only be used according to the user manual and the notes on use.

The information on setting values, application sites, duration of application and the use of the respective Erbe product is based on the clinical experience of physicians independent from Erbe. They represent guidelines which need to be checked by the user for their suitability for the actual planned application. Depending on the circumstances of an actual application case, it may be necessary to deviate from the information provided. The user has to check this on his/her own responsibility in each case when using an Erbe product. We wish to point out that science and technology is constantly subject to new developments arising from research and clinical experience. For this reason it may be necessary for the user to deviate from the information provided in this document.

This document contains information about Erbe products which may possibly not be approved in a specific country. The user of the respective Erbe product is obliged to inform him/herself whether the Erbe product he/she is using is legally approved in his/her country and/or if legal requirements or restrictions for use possibly exist and to which extent.

This document is not intended for users in the USA.

Erbe Elektromedizin GmbH
Waldhoernlestrasse 17
72072 Tuebingen
Germany

Phone +49 7071 755-0
Fax +49 7071 755-179
info@erbe-med.com
erbe-med.com