erbe
power your performance.

ENDO CUT® Q

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nível 3

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Literatura

- Ahmad NA, Kochman ML, Long WB, Furth EE, Ginsberg GG. Efficacy, safety, and clinical outcomes of endoscopic mucosal resection: a study of 1010 cases. *Gastrointestinal Endoscopy* 2002; 55(3):390-396.
- Binmoeller KF, Bohnacker S, Seifert H, Thonke F, Valdeyar H, Soehendra N. Endoscopic snare excision of "giant" colorectal polyps. *Gastrointest Endosc*. 1996; 43(3): 183-188.
- Brandimarte G. Endoscopic snare excision of large pedunculated colorectal polyps: A new, safe, and effective technique. *Endoscopy* 2001; 33(10):854-857.
- Brooker JC, Saunders BP, Shah SG, Williams CB. Endoscopic resection of large sessile colonic polyps by specialist and non-specialist endoscopists. *British Journal of Surgery* 2002; 89:1020-1024.
- Christie JP. Polypectomy or colectomy? Management of 106 consecutively encountered colorectal polyps. *Am Surg*. 1988;54(2):93-99.
- Classen M, Tytgat, GNJ, Lightdale CJ. *Gastroenterological Endoscopy*. Thieme-Verlag 2002.
- Deinlein P, Reulbach U, Stolte M, Vieth M, Risikofaktoren der lymphogenen Metastasierung von kolorektalen pT1-Karzinomen. *Pathologie* 2003;24:387-393.
- Di Giorgio P, De Luca L, Calcagno G, Rivellini G, Mandato M, De Luca B. Detachable snare versus epinephrine injection in the prevention of postpolypectomy bleeding: a randomized and controlled study. *Endoscopy* 2004; 36(10):860-863.
- Dobrowolski S, Dobosz M, Babicki A, Dymecki D, Hac S. Prophylactic submucosal saline-adrenaline injection in colonoscopic polypectomy: prospective randomized study. *Surgical Endoscopy* 2004; 18(6): 990-993.
- Dollhopf M, Heldwein W, Schmidtdorff G, Rösch T, Meining A, Schmitt W. Risk factors for perforation after polypectomy. Results of the Munich polypectomy study (MUPS) part 1. *Gastrointest Endosc* 2005; 55: AB99.
- Doniec JM, Lohnert MS, Schniewind B, Bokelmann F, Kremer B, Grimm H. Endoscopic removal of large colorectal polyps: prevention of unnecessary surgery? *Dis Colon Rectum*. 2003;46(3):340-348.
- Farin G. Möglichkeiten und Probleme der Standardisierung der Hochfrequenzleistung. In *Hochfrequenzdiathermie in der Endoskopie*; Springer-Verlag Berlin Heidelberg 1987.
- Fasoli R, Repaci G, Comin U, Minoli G, Italian Association of Hospital Gastroenterologists. A multi-center North Italian prospective survey on some quality parameters in lower gastrointestinal endoscopy. *Digestive and Liver Disease* 2002; 34(12):831-841.
- Frühmorgen P, Kriel L. Guidelines of the German Society of Digestive and Metabolic Diseases. Guidelines for endoscopic colorectal polypectomy with the sling. German Society of Digestive and Metabolic Diseases. *Z Gastroenterol*. 1998; 36(2):117-119.
- Kaneko E, Haradat H, Kasugi T, Sahita T. The result of a multicenter analysis from 1983-1992. *Gastroenterol. Endosc.* 1995;37:642-652.
- Garby JR, Suc B, Rotman N, Fourtainer G, Escat J. Multicentre study of surgical complications of colonoscopy. *British Journal of Surgery* 1996; 83:42-44.
- Giovannini M, Bories E, Presenti C, Moutardier V, Monges G, Danisi C, Lelong B, Delpero JR. Circumferential endoscopic mucosal resection in Barrett's esophagus with high-grade intraepithelial neoplasia or mucosal cancer. Preliminary results in 21 patients. *Endoscopy* 2004; 36(9):782-787.
- Giovannini M, Bernardini D, Moutardier V, Monges G, Houvenaeghel G, Seitz JF, Delpero JR. Endoscopic mucosal resection (EMR): results and prognostic factors in 21 patients. Preliminary results in 21 patients. *Endoscopy* 1999; 31(9):689-701.
- Greff M, Palazzo L, Ponchon TH, Canard JM; Council of The French Society of Gastrointestinal Endoscopy (S.F.E.D.). Guidelines of the French Society of Digestive Endoscopy: endoscopic mucosectomy. *Endoscopy* 2001; 33(2):187-190.
- Heldwein W, Dollhopf M, Schmidtdorff G, Rösch T, Meining A, Schmitt W. Risk factors for bleeding after polypectomy. Results of the Munich polypectomy study (MUPS) part 2. *Gastrointest Endosc* 2005; 55: AB99.
- Hurlstone DP, Sanders DS, Cross SS, Adam I, Shorthouse AJ, Brown S, Drew K, Lobo AJ. Colonoscopic resection of lateral spreading tumours: a prospective analysis of endoscopic mucosal resection. *Gut* 2004;53(9):1334-1339.
- Hsieh YH, Lin HJ, Tseng GY, Perng CL, Li AF, Chang FY, Lee SD. Is submucosal epinephrine injection necessary before polypectomy? A prospective, comparative study. *Hepatogastroenterology* 2001; 48(41): 1379-1382.
- Inoue H, Tani M, Nagai K, Kawano T, Takeshita K, Endo M, Iwai T. Treatment of esophageal and gastric tumors. *Endoscopy*. 1999; 31(1):47-55.
- Inoue H, Kawano T, Tani M, Takeshita K, Iwai T. Endoscopic mucosal resection using a cap: techniques for use and preventing perforation. *Can J Gastroenterol*. 1999; 13(6): 477-480.
- Kratzsch KH, Bosseckert H. Results of endoscopic polypectomy in the upper gastrointestinal tract-a multicenter study. *Deutsche Zeitschrift für Verdauungs- und Stoffwechselkrankheiten*.1984; 44(2):61-66.
- Löhnert MS, Wittmer A, Doniec JM. Endoscopic removal of large colorectal polyps. *Zentralbibliothek der Chirurgie* 2004; 129(4):291-295.
- Lux G, Semm K. *Hochfrequenzdiathermie in der Endoskopie*. Springer-Verlag Berlin Heidelberg, 1987 (Sonderdruck).
- Messmann H. *Lehratlas der Koloskopie*. Thieme Verlag 2004.
- Muehldorfer SM, Stolte M, Martus P, Hahn EG, Ell C. Diagnostic accuracy for forceps biopsy versus polypectomy for gastric polyps: a prospective multicenter study. *Gut* 2002; 50:465-470.
- Nivatvongs S. Complications in colonoscopic polypectomy. An experience with 1,555 polypectomies. *Dis Colon Rectum*. 1986;29(12): 825-830.
- Ono H, Kondo H, Gotoda T, Shirao K, Yamaguchi H, Saito D, Hosokawa K, Shimoda T, Yoshida S. Endoscopic mucosal resection for treatment of early gastric cancer. *Gut*. 2001; 48(2):225-229.
- Parra-Blanco A, Kaminaga N, Kojima T, Endo Y, Tajiri A, Fujita R. Colonoscopic polypectomy with cutting current: Is it safe? *Gastrointestinal Endoscopy* 2000; 51(6): 676-681.
- Parra-Blanco A, Kaminaga N, Kojima T, Endo Y, Uragami N, Okawa N, Hattori T, Takahashi H, Fujita R. Hemoclipping for postpolypectomy and postbiopsy colonic bleeding? *Gastrointestinal Endoscopy* 2000; 51(1):37-41.
- Ponchon T. Endoscopic mucosal resection. *J Clin Gastroenterol*. 2001; 32(1):6-10.
- Reed WP, Kilkenny JW, Dias CE, Wexner SD, SAGES EGD Outcomes Study Group. A prospective analysis of 3525 esophagogastroduodenoscopies performed by surgeons. *Surgical Endoscopy* 2004; 18:11-21.
- Rösch T, Sarbia M, Schumacher B, Deinert K, Frimberger E, Toerner T, Stolte M, Neuhaus H. Attempted endoscopic en bloc resection of mucosal and submucosal tumors using insulated-tip knives: a pilot series. *Endoscopy* 2004; 36(9):788-801.
- Schmitt W, Gospos J, Held T, Dollhopf M: Endoskopische Therapie der kolorektalen Neoplasie – neue Technik. *Dtsch med Wochenschr* 2000; 26:5250-5254.
- Schmitt W, Gospos J: Endoskopische Mukosektomie von breitbasigen großen Adenomen und T1-Karzinomen im Kolon. *Chir Gastroenterol* 2005; 21: 44-49
- Seitz U, Seewald S, Bohnacker S, Soehendra N. Advances in interventional gastrointestinal endoscopy in colon and rectum. *Int J Colorectal Dis*. 2003;18(1):12-18.
- Shim CS. Endoscopic mucosal resection: an overview of the value of different techniques. *Endoscopy*. 2001 Mar;33(3):271-275.
- Sieg A, Hachmoeller-Eisenbach U, Eisenbach T. Prospective evaluation of complications in outpatient GI endoscopy: a survey among German gastroenterologists. *Gastrointestinal Endoscopy* 2001; 53(6):620-627.
- Shim CS. Endoscopic mucosal resection: an overview of the value of different techniques. *Endoscopy* 2001; 33:271-275.
- Singh N, Harrison M, Rex DK. A survey of colonoscopic polypectomy practices among clinical gastroenterologists. *Gastrointestinal Endoscopy* 2004; 60(3):414-418.
- Takekoshi T, Baba Y, Ota H, Kato Y, Yanagisawa A, Takagi K, Noguchi Y. Endoscopic resection of early gastric carcinoma: results of a retrospective analysis of 308 cases. *Endoscopy* 1994; 26(4):352-358.
- Soehendra N, Binmoeller KF, Seifert H, Schreiber HW. *Praxis der therapeutischen Endoskopie*. Thieme Verlag 1997.
- Waye JD, Rex DK, Willima CB. *Colonoscopy*, 2003. Blackwell Publishing.
- Waye JD, Lewis BS, Frankel A, Geller SA. Small colon polyps. *Am J Gastroenterol* 1988; 83:120-122.
- Waye JD, Lewis BS, Yessayan S. Colonoscopy: a prospective report of complications. *Journal of Clinical Gastroenterology* 1992; 15(4):347-351.
- Waye JD, Kahn O, Auerbach ME. Complications of colonoscopy and flexible sigmoidoscopy. *Complications of Gastrointestinal Endoscopy* 1996; 6(2):343-377.

Erbe Elektromedizin GmbH
Waldhoernlestrasse 17
72072 Tuebingen
Alemanha

Tel +49 7071 755-0
Fax +49 7071 755-179
info@erbe-med.com
erbe-med.com