

No: 20195-230

BiClamp[®] 110



- Tissue structures can be coagulated over larger areas
- The individual dissection of vessels is usually not necessary
- Effective and rapid coagulation due to the optimal synchronization of the instrument and the VIO generator
- Highly economical: BiClamp is reusable; the procedure saves time and suture material
- AUTO STOP is user-friendly and offers optimal results
- When the instrument is plugged in, Plug & Play automatically adjusts the specific parameters

Dimensions and weight

Shaft diameter	5 mm
Length	110 mm
Cable length	with connecting cable 4 m

Basic features

Connection	MF plug
------------	---------

Surface

Branches	smooth
Coating	with non-adhesive coating

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Max. electrical capacity	0.25 kVp
Packaging unit	1 Piece