


No: 20402-037

Flexible Cryoprobe



Dimensions and weight

Dimensions	Ø 1,9 mm
Length	900 mm

Basic features

Flexibility	flexible
-------------	----------

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	Ila
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
----------------	---------