


No: 20150-015

HybridAPC



- Cost-efficient disposable instrument
- Multifunctional instrument with proven APC and waterjet functions
- One instrument for all Barrett localizations and esophagus lumina
- No instrument change required between elevation and ablation
- Easy to use
- Filter-integrated probe provides maximum protection against contamination

Dimensions and weight

Dimensions	ø 2,3 mm
Length	1.9 meter

Basic features

Compatibility	VIO, ERBEJET 2, APC 2
Flexibility	flexible
Connection	FiAPC plug, fluid connector

Cleaning and Reprocessing

Sterile packaging	EO
Single-use product	yes

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	5 Piece
Max. electrical capacity	5 kVp