

No: 20195-144

Insert for LAP forceps, fenestrated



Dimensions and weight

Needle length	340 mm
Shaft diameter	5 mm
Length	340 mm

Basic features

Shape / Forming	fenestrated
-----------------	-------------

Surface

Branches

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 134 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
Max. electrical capacity	0.25 kVp