

No: 20191-331



## Needle electrode

### Dimensions and weight

Shaft diameter	2 mm
Length	110 mm

### Basic features

Shape / Forming	needle angled 90°
-----------------	-------------------

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 134 °C
Single-use product	yes

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Max. electrical capacity	4 kVp
Packaging unit	5 Piece