

No: 20195-249

## BiClamp<sup>®</sup> E LAP forceps



- Effective and rapid coagulation due to the optimal synchronization of the instrument and the VIO generator
- Highly economical: BiClamp is reusable; the procedure saves time and suture material

### Dimensions and weight

Shaft diameter	5 mm
Length	340 mm
Cable length	with connecting cable 4 m

### Basic features

Shape / Forming	fenestrated
Connection	MF plug

### Surface

Branches	smooth
Coating	with non-adhesive coating

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	1 Piece
Max. electrical capacity	0.25 kVp