

No: 21191-361

## Spatula electrode



### Dimensions and weight

Dimensions	2,3 x 19 mm
Shaft diameter	2.4 mm
Length	45 mm

### Basic features

Shape / Forming	straight
Colour	blue

### Surface

Insulation	insulated
------------	-----------

### Cleaning and Reprocessing

Sterile packaging	E0
Single-use product	yes

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	20 Piece
Max. electrical capacity	4.3 kVp