

No: 20195-554

## Bipolar forceps Premium, bayonet

- Limited tissue adhesion due to high temperature conductivity
- Long instrument life due to high-quality alloy
- Targeted precise coagulation due to precise design of the tips
- Firm grasping of tissue due to abraded gripping surface
- Less time required intraoperatively for cleaning
- Good visibility at the operating site due to reduced light reflections (e.g. during preparation under a microscope)



### Dimensions and weight

Tip	tip 1 mm, blunt
Length	250 mm

### Basic features

Shape / Forming	bayonet
-----------------	---------

### Surface

Branches	roughened
----------	-----------

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	1 Piece
Max. electrical capacity	0.5 kVp