

No: 10140-500

VIO[®] 100 C



- Intuitive operation
- Maximum CUT output 100 watt
- Configurable sockets (monopolar, bipolar, neutral)
- AUTO START
- 4 Predefined programs (overwriteable)
- Applications: interventions performed on an outpatient basis in the fields of dermatology, plastic surgery, ENT, gynecology, casualty surgery, ophthalmology

Power output

Maximum CUT output	up to 100 watt
Maximum COAG output	up to 80 watt
Safety system	yes
Frequency	510-550 kHz

Mains connection

System voltage	100-120 / 220-240 Volt
Power frequency	50 / 60 Hz
Power input during stand-by	15 Watts
Power input during max RF output	[190 Watt / 200 VA]
Potential equalization connection	yes

Dimensions and weight

Width x Height x Depth	280 x 135 x 300 mm
Weight	4 kg

Additional product information

Packaging unit	1 Piece
----------------	---------

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF
Category acc. to UMDNS	