

No: 20195-035



## Bipolar forceps Classic, bent

### Dimensions and weight

Tip	tip 0.5 mm, fine
Length	105 mm
Gripping length	6 mm

### Basic features

Shape / Forming	bent
-----------------	------

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Max. electrical capacity	0.3 kVp
Packaging unit	1 Piece