

No: 21191-312

Needle electrode



Dimensions and weight

| | |
|----------------|---------------|
| Dimensions | ø 0,8 x 22 mm |
| Shaft diameter | 2.4 mm |
| Length | 40 mm |

Basic features

| | |
|-----------------|----------|
| Shape / Forming | straight |
| Colour | blue |

Cleaning and Reprocessing

| | |
|--------------------|-----|
| Sterile packaging | EO |
| Single-use product | yes |

Standards

| | |
|---|-----|
| Classification acc. to 93/42 EWG | IIb |
| Protection category acc. to EN 60 601-1 | |
| Type acc. to EN 60 601-1 | |

Additional product information

| | |
|--------------------------|----------|
| Packaging unit | 50 Piece |
| Max. electrical capacity | 4.3 kVp |