

No: 21191-312

## Needle electrode



### Dimensions and weight

Dimensions	ø 0,8 x 22 mm
Shaft diameter	2.4 mm
Length	40 mm

### Basic features

Shape / Forming	straight
Colour	blue

### Cleaning and Reprocessing

Sterile packaging	EO
Single-use product	yes

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	50 Piece
Max. electrical capacity	4.3 kVp