

No: 21191-468

Coated Needle Electrode



Dimensions and weight

Dimensions	ø 0,9 x 16 mm
Shaft diameter	2.4 mm
Length	50 mm

Basic features

Shape / Forming	straight
Colour	blue

Surface

Insulation	insulated
Coating	[PTFE-Beschichtung]

Cleaning and Reprocessing

Sterile packaging	EO
Single-use product	yes

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	10 Piece
Max. electrical capacity	4 kVp