

No: 21191-457

Coated Spatula Electrode



Dimensions and weight

Dimensions	2,3 x 19 mm
Shaft diameter	2.4 mm
Length	134 mm

Basic features

Shape / Forming	straight
Colour	blue

Surface

Coating	[PTFE-Beschichtung]
---------	---------------------

Cleaning and Reprocessing

Sterile packaging	EO
Single-use product	yes

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	10 Piece
Max. electrical capacity	4 kVp