

No: 21191-013

Needle electrode

- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines



Dimensions and weight

Dimensions	∅ 0,8 x 29 mm
Shaft diameter	4 mm
Length	65 mm

Basic features

Shape / Forming	straight
Compatibility	VIO, ICC, ACC
Flexibility	rigid
Colour	black

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
Max. electrical capacity	4.3 kVp