

No: 20190-170

Electrosurgical pencil with rocker switch



Dimensions and weight

| | |
|----------------|---------------------------|
| Shaft diameter | 2.4 mm |
| Cable length | with connecting cable 5 m |

Basic features

| | |
|---------------|-------------------------------|
| Compatibility | VIO, ICC, ACC, non-ERBE units |
| Connection | International |

Cleaning and Reprocessing

| | |
|--------------------|-----|
| Sterile packaging | E0 |
| Single-use product | yes |

Standards

| | |
|---|-----|
| Classification acc. to 93/42 EWG | IIb |
| Protection category acc. to EN 60 601-1 | |
| Type acc. to EN 60 601-1 | |

Additional product information

| | |
|--------------------------|----------|
| Packaging unit | 25 Piece |
| Max. electrical capacity | 5 kVp |