


No: 20402-410

## Flexible cryoprobe, single-use



### Dimensions and weight

Dimensions	ø 1,7 mm
Length	1150 mm

### Basic features

Flexibility	flexible
-------------	----------

### Cleaning and Reprocessing

Sterile packaging	E0
Single-use product	yes

### Standards

Classification acc. to 93/42 EWG	Ila
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	5 Piece
----------------	---------