

No: 10321-000

IES 300



- Activation either manually, by foot switch or automatically by an activation sensor coupled to the HF-generator
- High flow rate yet quiet operation
- Easy-to-use control panel; easy wipe disinfection allows maximum hygiene
- Comprehensive safety concept: e.g. automatic evacuation stop in the event of tissue suction (during high flow rate)
- Cost-efficient: long filter life (35 h) under max. flow rate

Mains connection

Power frequency	50 / 60 Hz
-----------------	------------

Dimensions and weight

Width x Height x Depth	410 x 210 x 368 mm
Weight	14 kg

Additional product information

Packaging unit	1 Piece
----------------	---------

Standards

Classification acc. to 93/42 EWG	I
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF