

No: 10140-400

VIO[®] 200 S



- Automatic output dosage for all regulative technologies
- Precise power output adjusted to the individual indication
- Variable configuration with socket modules
- Simple to operate due to FocusView
- Display-supported neutral electrode safety system NESSY provides interactive assistance when positioning the return electrode
- Up to 10 programs can be stored
- 2 footswitches can be connected to the unit
- VIO 200 S can be updated (updates possible)

Power output

Maximum CUT output	200 watt at 500 ohm
Maximum COAG output	up to 120 watt
Safety system	NESSY
Frequency	350 kHz

Mains connection

Power frequency	50 / 60 Hz
Power input during stand-by	40 Watts
Power input during max RF output	500 watt / 920 VA
System voltage	100-120 / 220-240 Volt
Potential equalization connection	yes

Dimensions and weight

Width x Height x Depth	410 X 165 X 380 mm
Weight	8,8 kg

Additional product information

Packaging unit	1 Piece
----------------	---------

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF
Category acc. to UMDNS	