

No: 10140-200

VIO<sup>®</sup> 200 D



- Automatic output dosage for all regulative technologies
- Precise power output adjusted to the individual indication
- New and improved CUT and COAG functions
- Variable configuration with socket modules
- Upgrades for specific medical specialties are possible at any time
- Simple, interactive and safe operation using the TFT display and FocusView
- Display-supported neutral electrode safety system NESSY provides interactive assistance when positioning the return electrode
- Individually configurable programs in simple text: switch on the unit and begin working
- ReMode function: "remote control" via the handle or foot switch straight from the sterile field
- Designed to provide solutions for specific medical indications in the medical fields gastroenterology, ENT, and bronchoscopy, and for use in operative outpatient procedures

#### Power output

Maximum CUT output	200 watt at 500 ohm
Maximum COAG output	up to 120 watt
Safety system	NESSY
Frequency	350 kHz

#### Mains connection

Power frequency	50 / 60 Hz
Power input during stand-by	40 Watts
Power input during max RF output	500 watt / 920 VA
System voltage	100-120 / 220-240 Volt
Potential equalization connection	yes

#### Dimensions and weight

Width x Height x Depth	410 x 165 x 380 mm
Weight	9,5 kg

#### Additional product information

Packaging unit	1 Piece
----------------	---------

#### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF
Category acc. to UMDNS	11-490