

No: 20196-119

Bipolar connecting cable



Dimensions and weight

Cable length	[4 m]
--------------	-------

Basic features

Compatibility	for WOLF resectoscopes
Connection	[mit MF-Stecker]

Cleaning and Reprocessing

Sterilization in autoclave	max 134 °C
Cleaning in washer disinfectant	max 95 °C

Standards

Classification acc. to 93/42 EWG	I
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Max. electrical capacity	1 kVp
Packaging unit	1 Piece