

No: 20196-118

Bipolar connecting cable



Dimensions and weight

| | |
|--------------|-------|
| Cable length | [4 m] |
|--------------|-------|

Basic features

| | |
|---------------|-------------------------|
| Compatibility | for STORZ resectoscopes |
| Connection | [mit MF-Stecker] |

Cleaning and Reprocessing

| | |
|---------------------------------|------------|
| Sterilization in autoclave | max 134 °C |
| Cleaning in washer disinfectant | max 95 °C |

Additional product information

| | |
|--------------------------|---------|
| Max. electrical capacity | 1 kVp |
| Packaging unit | 1 Piece |