

No: 20194-077

## Return electrode cable



### Dimensions and weight

Cable length	[4 m]
--------------	-------

### Basic features

Compatibility	VIO, ICC, ACC, T-Series
---------------	-------------------------

Connection	Standard
------------	----------

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
---------------------------------	-----------

Sterilization in autoclave	max 138 °C
----------------------------	------------

### Standards

Classification acc. to 93/42 EWG	I
----------------------------------	---

Protection category acc. to EN 60 601-1	
---	--

Type acc. to EN 60 601-1	
--------------------------	--

### Additional product information

Packaging unit	1 Piece
----------------	---------

Max. electrical capacity	0.5 kVp
--------------------------	---------